

# Performer Dental Chair and Modules Installation Checklist

For more detailed information, see the installation guides for the dental chair and optional modules.

	Model	S/N	Model	S/N
Chair:	Performer		Light:	
Delivery:			Assistant's:	Performer
Support Center:	Performer		Other:	

## Dental Chair

- Dental chair is securely anchored to the floor
- Chair functions properly, including the chair's lift, tilt, and programmable functions controlled by the touchpads and footswitch
- Headrest functions properly
- Stop plate functions properly

## Delivery System

- Handpieces function properly and are set to doctor preferences:
  - Water coolant, air coolant, and spray
  - Handpiece tubing flush
  - Handpiece holder valves and switches
  - Handpiece drive air pressure
- Performer Delivery System only: Rotation stop pin is installed and functioning properly
- Performer Delivery System only: Height of delivery system and rotation tension are properly adjusted
- 332/333 Delivery System only: Flexarm air brake functions properly
- 332/333 Delivery System only: Spring-assisted flexarm is counterbalanced properly
- Control head is level as measured on tray
- Control head rotation tension is properly adjusted
- Articulating arms do not drift
- Foot control operates properly

## Dental Light

- All light settings and controls work properly
- Dental light flexarm and head tension are properly adjusted

## Assistant's Instrumentation, Cuspidor, and Cupfill

- Vacuum/suction is set to the following specifications:
  - Wet vacuum** -  $10 \pm 2$  inches of Hg ( $34 \pm 7$  kPa), 9 scfm (255 SL/min) minimum
  - Dry/semi-dry vacuum** -  $4.5 \pm 1$  inches of Hg ( $16 \pm 3.5$  kPa), 12 scfm (340 SL/min) minimum
- Air and water syringe buttons operate smoothly
- Valves on HVE and saliva ejector move freely
- Cupfill timing is set properly
- Cuspidor flow pattern effectively rinses the bowl and bowl drains properly

## Utilities

- Regulated air pressure is at 80 psi (552 kPa)
- Floor box is free of air and water leaks
- The gravity drain functions properly
- Tubing and connections are not kinked and are free of air and water leaks
- Excess tubing and wires are coiled and stored away from moving parts

---

**Doctor / Clinic:**

---

**Address:**

---

**Telephone:**

---

**Installation Date:**

---

**Dealer / Branch:**

---

**Installer:**

---

**Operatory:**

---

**Notes:**

---

**Reviewed with doctor / customer:**

\_\_\_\_\_  
Signature

\_\_\_\_\_  
Date

**Doctor / customer unavailable to review checklist**



2601 Crestview Drive  
Newberg, OR 97132 USA  
Tel: 1.800.547.1883 Within USA/CAN  
Tel: +1.503.538.7478 Outside USA/CAN  
www.a-dec.com

*A-dec Inc. makes no warranty of any kind with regard to the content in this document including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose.*



86.0657.00 Rev C  
Date of Issue 2020-03-02  
Copyright 2020 A-dec Inc.  
All rights reserved.  
IGrphpor11