



With ICX, you have
peace of mind with your
waterline maintenance.



A-dec® ICX

Easy waterline maintenance. No mixing or measuring.

A-dec ICX waterline treatment tablets offer a unique and effective way to maintain clean dental unit waterlines—and it's easy to use since the effervescent ICX tablet does the work for you.* ICX continues working in the system with each treatment, protecting the waterlines from further contamination and reducing the need for frequent shock treatments.

- Simply add one tablet to an empty self-contained water bottle before each filling.
- As the tablet dissolves, it releases ingredients into the water that help maintain less than or equal to 10 colony forming units (CFU) per milliliter.
- ICX prevents the accumulation of odor and foul-tasting bacteria.
- Convenient single-dispensing packets eliminate measuring, mixing and messy solutions. Available for 0.7-liter and 2-liter bottles

To learn more about the A-dec ICX waterline maintenance, contact your local authorized dealer.

A-dec ICX

With ICX, your waterline maintenance protocol is not only easy and effective, but also safe for your investment since the concentration of ingredients in the water remains very low.

Convenient

- Specially formulated effervescent tablets
- Simple single-dose package; drop tablet into an empty self-contained water bottle and fill

Effective residual protection

- Maintains dental unit waterline effluent ≤ 10 CFU/ml
- Eliminates odor and foul-tasting bacteria
- Reduces build-up of deposits
- Continuously present in water system; provides a preventive, proactive (vs. reactive) solution
- Remains in water up to two weeks during periods of equipment non-use
- Eliminates the need to purge waterlines at night
- EPA registered
- Peer-reviewed studies on ICX:*
 - Journal of Clinical Dentistry JCD 2004, Vol. XV, Number 1, pp. 17-21
 - Journal of the American Dental Association, JADA June 2004, Vol. 135, pp. 799-805

Safety**

- Non-corrosive to dental equipment
- No odor or taste
- Includes three primary active ingredients: sodium percarbonate, silver nitrate, and cationic surfactants

Distribution, packaging and use

- Distributed exclusively through A-dec's worldwide authorized dealer network
- Tablets are available for two bottle sizes: 0.7- and 2-liter
- Individually foil-wrapped tablet, packaged 50 per box
- Cost effective; low annual cost
- Approximate usage rates for a single treatment room (one 2-liter or one and a half 0.7-liter bottles per day for a 16-day month):
 - Two months per 50-count box, 0.7-liter tablet
 - Three months per 50-count box, 2-liter tablet

Bonding

- Compatible for use with restorative dental materials; does not affect bond strength
- Clinical study on dentin bonding:*
 - Journal of Clinical Dentistry, JCD 2004, Vol. XV, Number 1, pp. 28-32

Product Description	Bottle Size	Contents	Total Treatments	A-dec Part Number
ICX Tablets	2-liter	1 case containing 36 boxes: 50 tablets per box	1,800 treatments	90.1065.00
ICX Tablets	0.7-liter	1 case containing 36 boxes: 50 tablets per box	1,800 treatments	90.1064.00

To learn more about the A-dec ICX waterline maintenance, contact your local authorized dealer. Visit a-dec.com/ICX or ask your authorized dealer for more information.

* For independent clinical study reprints or additional information, contact A-dec Customer Service at 1.800.547.1883 (1.503.538.7478 worldwide).

** Safety Data Sheet (SDS) are available at a-dec.com.

A-dec, the A-dec logo, A-dec ICX, are trademarks of A-dec Inc. and are registered in the United States and other countries. None of the trademarks or trade names in this document may be reproduced, copied, or manipulated in any manner without the express, written approval of the trademark owner.

© 2022 A-dec Inc. All rights reserved. 85.6093.00/1M/9-22/RevL