

A-dec Headquarters
 2601 Crestview Drive Newberg, Oregon 97132 USA
 Tel: 1.800.547.1883 Within USA /Canada
 Tel: +1.503.538.7478 Outside USA /Canada
www.a-dec.com

**A-dec 311ft/311B
 Dental Chair
 Mini Template**

86.0551.00 Rev E Copyright 2023 A-dec, Inc.
All rights reserved.

Date of Issue 2023-11-16

Plumbing Elevations:



1" PVC (1.315" O.D.) Vacuum

1-1/2" (38 mm) nominal pipe for gravity drain (when required)

3-3/4" (95 mm) max. height

1" to 2" (25 - 50 mm)

1" (25 mm) finish stub height

5" (127 mm) max. height

1/2" (13 mm) nominal pipe compression fitting with right-angle manual shut-off valves with 3/8" compression outlet for air and water (supplied by contractor)

PLEASE NOTE: **This drawing is for reference purposes ONLY. It is NOT TO SCALE.**